



AEHia



AEHIS



AEHIT



FALL
SUMMIT
2019

Building a Connected Asset Security Program: Lessons From the Field

Rodney Graves, GBMC

Toby Gouker, Carter Groome, First Healthcare
Advisory Solutions

FALL
SUMMIT
2019

About GBMC...

- Established in 1965 to combine the best of community and university-level medicine
- 342 bed medical center that that handles more than 23,000 admissions and over 52,000 emergency room visits annually
- GBMC Health Partners in Hunt Manor, Hunt Valley, Owings Mills, Texas Station, Perry Hall, Joppa Road and Jarrettsville
- Gilchrist hospice serves more than 7,600 patients each year



FALL
SUMMIT
2019

About First...

- Asset Orchestration & Risk Management
- Full Suite Privacy & Security Services
- Clinical and Revenue Cycle HIT



FALL
SUMMIT
2019

About GBMC Cybersecurity...

- 2016 Medstar, 2017 WannaCry & Not Petya
- Officially began with dedicated staff in 2016
- NIST Framework
- Emerging threat: IOT



FALL
SUMMIT
2019

Search for a Solution

- Full network visibility
- Internet research, industry colleagues and vendors
- Partnership Healthcare focused IT security



Program, not Product Solution



Technical and Organizational Change Solution

*Connected asset security is a classic technical risk management issue
With tricky organizational dynamics issues*

RISK MANAGEMENT

- Assess
- Detect
- Respond
- Defend/Protect



ORGANIZATIONAL DYNAMICS


- Attachment to current standard ways/culture
- Resistance to change
- Lack of IoT security education

01 LAYER

GAIN VISIBILITY
Inventory, calculate risk, and prioritize



- Develop relationship between IT, Security and HTM teams
- Coordinated active sample scan
- Discovery of any significant issues
- Initiate system wide passive scanning
- Conduct a risk-based assessment
- Set critical alert mitigation priorities



FALL
SUMMIT
2019


01 LAYER

GAIN VISIBILITY
Inventory, calculate risk, and prioritize



Platform initiation

- Manage installation, testing and operation
- Inventory connected assets within assigned network
- Develop a prioritized list of issues and vulnerabilities
- Prioritize for FY 2020 short-term plan and budget to address mission critical remediations



FALL
SUMMIT
2019

02 LAYER

OPERATIONALIZE Data-driven action plan to mitigate risk



- Inventory and Vulnerability verification
- Dashboard and report design and customization
- Mitigation prioritization and planning
- Roadmap and budget plan for connected asset program

Impact Risk	Problem / Vulnerability	Action Steps / Mitigation Choice	Status of Action / Next Steps
Critical	25MASTER Nurse Call System - SMB v1 Use	Microsoft Security Updates Needed	Network segmentation for Nurse Call Systems
Critical	27MASTER Nurse Call System - SMB v1 Use	Microsoft Security Updates Needed	Network segmentation for Nurse Call Systems
Critical	36MASTER Nurse Call System - SMB v1 Use	Microsoft Security Updates Needed	Network segmentation for Nurse Call Systems
Critical	46MASTER Nurse Call System - Windows XP	Replace EOL/EOS Microsoft XP Systems	Network segmentation for Nurse Call Systems
Critical	Axis IP Cameras - Access Control Bypass	Apply latest firmware release; Apply ACL to Axis Cameras	owned by SIM lab - Reached out
Critical	GBMCSEND JBoss JMXInvokerServlet remote code execution attempt	Apply latest security updates	Device located on RADNET. Working with personnel to identify
Critical	Windows Computers - RDP vulnerability	Apply latest security updates	patch upgrade schedule - security?
Critical	RDP vulnerability - APP18 & APP19	Servers identified for replacement/retirement	Timeline? Other compensating controls?
Critical	Windows Computers - SMB v1 vulnerabilities	Apply latest security updates	patch upgrade schedule - security?
Critical	SAFETYNET1 Patient Monitoring - SMB v1 Use	Microsoft Security Updates Needed	patch upgrade schedule - security?
Critical	SAFETYNET2 Patient Monitoring - SMB v1 Use	Microsoft Security Updates Needed	patch upgrade schedule - security?
High	77 Hospira Infusion pumps - DoS vulnerability	Apply latest firmware release	emailed HTM to get confirmation
High	77 Hospira Infusion pumps - No root telnet authentication required	Apply latest firmware release	emailed HTM to get confirmation
High	77 Hospira Infusion pumps - No network data validation	Apply latest firmware release	emailed HTM to get confirmation
High	77 Hospira Infusion pumps - Remote attacks via code execution	Apply latest firmware release	emailed HTM to get confirmation
High	SIEMENS F2 - FTP login using compromised password	Update login credentials and use secure FTP	Has this been completed? (GBMC)
High	Zebra Label Printers using default login credentials	Update login credentials	Vendor is providing pricing for a management tool to identify and change credentials.
High	MININT-3GBOEA4 - SNMP v1 communication with unknown devices	Validate SNMP requirements and update to SNMP v3, if applicable	Scheduled for upgrade on 9/10/2019 per Barry Bogardus email on 9/8/2019

02 LAYER

OPERATIONALIZE
Data-driven action plan to mitigate risk



- Mitigation Program Activities
 - Network segmentation
 - Subnet reassignments
 - Patches
 - Equipment retirement
 - System hardening
 - Compensating Controls...

Submission Information			
Date		Submitted By	
Vendor Name		Vendor Category	
+			
Device and Technical Identification Details			
Device Name		Device Domain Name	
Device IP Address		Device IP Subnet	
Device MAC Address		Device OS	
OS Service Pack Level			
Security and Vulnerability Details and Information			
Security Alert Type		Security CVE ID	
Impact Risk		Remotely Exploitable	
Security Vulnerability Description			
Remediation Action Plan and Related Details			
Remediation Timeline		Change Submittal Date	
Change Approval Date		Change Tracking Number	
Implement Date		Implementation Result	
Remediation Plan Details			
Post Implementation Actions and Notes			
Post Action Security Scan Date		Scan Results	
Final Closure Date		Signature	
Post Closure Final Notes			

FALL
SUMMIT
2019

03
LAYER

GOVERN

Security vigilance through policy,
people, and lifecycle management




Communication to C-suite
Policies for new purchases
Committee to handle change
management
Monitor regulatory changes &
manufacturer alerts
Business operation changes
Closed-loop equipment lifecycle
management
Continuous monitoring

Organizational Cyber Health Maturity Level Attainment Objectives	Weighted %
Visibility, inventory, workforce education and awareness	15%
24x7 monitoring and alerting	10%
Full system integration (CMMS, SIEM, NAC, vulnerability scanners, firewalls, SOAR, Network management, wireless LAN, managed security services)	10%
Complete mitigation of critical and high vulnerabilities	25%
Mitigation of remaining issues to attain risk appetite	15%
Fully established governance program – <ul style="list-style-type: none">• Data driven operationalization of policies and processes• Enhanced inventory capabilities for utilization metrics, device load leveling, and maintenance	25%
Total	100%

FALL
SUMMIT
2019

03 LAYER

GOVERN Security vigilance through policy, people, and lifecycle management



03
LAYER

GOVERN
Security vigilance through policy,
people, and lifecycle management



Network Usage

ESTABLISHED CONNECTIONS

Continuous Monitoring


Chart Type: Category View Detail View

X-Ray GE DR6000

Geo-location

IP/Users

Applications




FALL
SUMMIT
2019



Continuous Monitoring Bonuses

- Capacity utilization metrics for new purchase justification
- Device load leveling to avoid new purchases
- Maintenance based on use, not time passed
- Abnormalities as indication of pending failures
- Performance improvement coaching
- Redirect of patient to available capacity



Questions & Answers

Now, or later....

rgraves@gbmc.org

tgouker@fcp.com

cgroome@fcp.com